National Diagnostic Protocol

Phyllosticta ampelicida

the cause of black rot on grapevine



NDP 13 V2

© Commonwealth of Australia

Ownership of intellectual property rights

Unless otherwise noted, copyright (and any other intellectual property rights, if any) in this publication is owned by the Commonwealth of Australia (referred to as the Commonwealth).

Creative Commons licence

All material in this publication is licensed under a Creative Commons Attribution 3.0 Australia Licence, save for content supplied by third parties, logos and the Commonwealth Coat of Arms.



Creative Commons Attribution 3.0 Australia Licence is a standard form licence agreement that allows you to copy, distribute, transmit and adapt this publication provided you attribute the work. A summary of the licence terms is available from http://creativecommons.org/licenses/by/3.0/au/deed.en. The full licence terms are available from http://creativecommons.org/licenses/by/3.0/au/deed.en. The full licence terms are available from https://creativecommons.org/licenses/by/3.0/au/deed.en.

This publication (and any material sourced from it) should be attributed as: Subcommittee on Plant Health Diagnostics (2017). National Diagnostic Protocol for *Phyllosticta ampelicida* – NDP13 V2. (Eds. Subcommittee on Plant Health Diagnostics) Authors Sosnowski, M. and Tan, Y-P; Reviewers de Alwis, S.K. Tan, Y-P. and Baskarathevan, J. ISBN 978-0-9945113-7-9. CC BY 3.0.

Cataloguing data

Subcommittee on Plant Health Diagnostics (2017). National Diagnostic Protocol for *Phyllosticta ampelicida* – NDP13 V2. (Eds. Subcommittee on Plant Health Diagnostics) Author Sosnowski, M and Tan, Y-P; Reviewers de Alwis, S.K., Tan, Y-P. and Baskarathevan, J. ISBN 978-0-9945113-7-9

ISBN 978-0-9945113-7-9

Internet

Report title is available at: <u>http://plantbiosecuritydiagnostics.net.au/resource-hub/priority-pest-diagnostic-resources/</u>

Department of Agriculture and Water Resources

Street Address 1: <u>18 Marcus Clarke Street</u>, Canberra City ACT 2601
Street Address 2: <u>7 London Circuit</u>, Canberra City ACT 2601
Postal Address: GPO Box 858, Canberra City ACT 2601
Switchboard Phone: 02 6272 3933
Web: <u>http://www.agriculture.gov.au</u>

Inquiries regarding the licence and any use of this document should be sent to: <u>copyright@agriculture.gov.au</u>.

The Australian Government acting through the Department of Agriculture and Water Resources has exercised due care and skill in the preparation and compilation of the information and data in this publication. Notwithstanding, the Department of Agriculture and Water Resources, its employees and advisers disclaim all liability, including liability for negligence, for any loss, damage, injury, expense or cost incurred by any person as a result of accessing, using or relying upon any of the information or data in this publication to the maximum extent permitted by law.

Purpose

National Diagnostic Protocols (NDPs) are diagnostic protocols for the unambiguous taxonomic identification of plant pests. NDPs:

- are a verified information resource for plant health diagnosticians
- are consistent with ISPM No. 27 Diagnostic Protocols for Regulated Pests
- provide a nationally consistent approach to the identification of plant pests enabling transparency when comparing diagnostic results between laboratories; and,
- are endorsed by regulatory jurisdictions for use (either within their own facilities or when commissioning from others) in a pest incursion.

Where an International Plant Protection Convention (IPPC) diagnostic protocol exists it should be used in preference to NDPs although NDPs may contain additional information to aid diagnosis. IPPC protocols are available on the IPPC website:

https://www.ippc.int/core-activities/standards-setting/ispms

Process

NDPs are facilitated and endorsed by the Subcommittee on Plant Health Diagnostics (SPHD). SPHD reports to Plant Health Committee and is Australia's peak technical and policy forum for plant health diagnostics.

NDPs are developed and endorsed according to Reference Standards developed and maintained by SPHD. Current Reference Standards are available at <u>http://plantbiosecuritydiagnostics.net.au/sphd/sphd-reference-standards/</u>

NDPs are living documents. They are updated every 5 years or before this time if required (i.e. when new techniques become available).

Document status

This version of the National Diagnostic Protocol (NDP) for *Phyllosticta ampelicida* is current as at the date contained in the version control box below.

PEST STATUS	Not present in Australia
PROTOCOL NUMBER	NDP 13
VERSION NUMBER	V2
PROTOCOL STATUS	Endorsed
ISSUE DATE	2017
REVIEW DATE	2022
ISSUED BY	SPHD

The most current version of this document is available from the SPHD website: http://plantbiosecuritydiagnostics.net.au/resource-hub/priority-pest-diagnostic-resources/

Further information

Inquiries regarding technical matters relating to this project should be sent to: <u>sphds@agriculture.gov.au</u>

Contents

1	INTRO	DDUCTION
	1.1	Host range
2	TAXO	NOMIC INFORMATION
	2.1	Taxonomic description
	2.2	Morphological description4
3	DETE	CTION
	3.1	Leaf, stem and fruit symptoms6
	3.2	Confusion with other diseases10
	3.3	Sampling methods
4	IDENT	TIFICATION12
	4.1	Isolation of pathogen from symptomatic plant material12
	4.2	Morphological examination
	4.3	Molecular method16
5	CONT	ACTS FOR FURTHER INFORMATION19
6		OWLEDGEMENTS
7	REFE	RENCES
	7.1	Other useful references
8	APPE	NDICES
	8.1	Disease cycle

1 INTRODUCTION

Black rot of grape is caused by the ascomycete fungus *Phyllosticta ampelicida* (syn. *Guignardia bidwellii*) and can lead to substantial yield losses in humid regions when it is not effectively managed. Most cultivars of *Vitis vinifera*, as well as French/American hybrids and American bunch grapes, are susceptible, while varieties of muscadine range in disease susceptibility from resistant to very susceptible.

The fungus may be spread on grapevine cuttings and fruit.

1.1 Host range

Sivanesan A and Holliday P (1981)

1.1.1 Primary hosts

Vitis arizonica	Canyon grape
Vitis labrusca	American grape
Vitis rotundifolia	Muscadine grape
Vitis vinifera	Domestic grape

1.1.2 Alternative hosts

Wild grape
Birds nest fern
Ornamental vine
Virginia creeper
Boston ivy
Amur grape

2 TAXONOMIC INFORMATION

2.1 Taxonomic description

Class: Dothideomycetes *Order* Incertae sedis *Family*: Phyllostictaceae *Genus*: Phyllosticta *Species*: ampelicida

Name

Phyllosticta ampelicida (Engelm.) Aa, Stud. Mycol. 5: 28 (1973)

Synonyms

Botryosphaeria bidwellii (Ellis) Petr. 1958 Carlia bidwellii (Ellis) Magnus 1892 Carlia bidwellii (Ellis) Prunet 1898 Depazea labruscae Engelm. 1877 Guignardia bidwellii (Ellis) Viala & Ravaz 1892 Guignardia bidwellii f.sp. bidwellii (Ellis) Viala & Ravaz 1892 Guignardia bidwellii var. bidwellii (Ellis) Viala & Ravaz Guignardia bidwellii var. euvitis Luttr. 1946 Guignardia bidwellii f. parthenocissi Luttr. 1946 Laestadia bidwelli (Ellis) Viala & Ravaz, 1888 Naemospora ampelicida Engelm. 1863 Phoma ustulata Berk. & M.A. Curtis 1873 Phoma uvicola var. labruscae Thüm. 1878 Phyllachorella bidwellii (Ellis) Theiss. 1919 Phyllosticta ampelopsidis Ellis & G. Martin 1886 Phyllostictina clemensae Petr. 1928 *Phyllosticta labruscae* Thüm. 1878 Phyllosticta labruscae var labruscae Thüm. 1878 Phyllosticta vitea Sacc. 1896 Phyllosticta viticola Thüm. 1878 Physalospora bidwelli (Ellis) Sacc., 1882 Sphaeralla bidwellii (Ellis) Ellis 1890 Sphaeria bidwelli Ellis 1880 Septoria viticola

2.2 Morphological description

The following description is taken from Sivanesan and Holliday (1981).

Pseudothecia formed as locules in a stroma, depressed globose, immersed subepidermal, 70-180 μ m broad with a flat or papillate ostiolar apex (Fig. 1A). The pseudothecial wall is made up of pseudoparenchymatic cells. Asci arise from a cushion shaped hyaline tissue at the base, cylindrical to clavate, 45-65 × 9-14 μ m (Fig. 1B). Ascospores hyaline, one celled, ovoid to ellipsoid, 12-17 × 6-7.5 μ m, often with hyaline, mucilaginous, apical caps.

Pycnidia mostly epiphyllous, solitary, unilocular, globose or depressed globose, with a flat or inconspicuous papillate ostiolar apex, 120-230 μ m broad (Fig. 1C). Stroma poorly developed on leaves but well developed on fruits. Conidiogenous cells conical to cylindrical (Fig 1D). Conidia one celled, hyaline, broadly ovoid, ellipsoidal or almost globose, somewhat clavate when young and slightly indented, 5-12 × 4-7 μ m, surrounded by a mucilaginous sheath and with an apical hyaline appendage as long as the conidium. Spermatia hyaline, unicellular, rod shaped 4-7 × 0.5-2 μ m.




Figure 1. Microscopic description of *Guignardia bidwellii*; A. pseudothecium, B. ascus and ascospores, C. pycnidium and D. conidiogenous cells and conidia (Sivanesan and Holliday 1981).

Descriptions from Zhang et al. (2013) based on neotype.

Pycnidia black or blackish brown, globose or subglobose, 150-300 mm diam. on PDA media. Conidiogenous cells cylindrical or conical, $9-13(-15) \times 3-4.5 \text{ mm}$. Conidia hyaline one-celled, broadly ovoid, ellipsoidal, seldomsubglobose, some with a truncate base, $8.5-12.5 \times 6-7.5 \text{ mm}$, surrounded by a slime layer, containing small green guttules, with short appendages, 2-7 mm. Spermatia hyaline, one-celled, dumbbell shaped, $5-7 \times 1.3-2.2 \text{ mm}$.

Culture characteristics: Colonies reaching 3 cm diam. in 7 d on PDA, pale gray to olive green when mature, gray olivaceous to greenish black in reverse. Aerial mycelium dense, white; submerged mycelium olivaceous to brown. Stromata formed after 20 d, black.

3 DETECTION

Phyllosticta ampelicida can be readily identified by macroscopic symptoms on most parts of the plant. The disease cycle of black rot is described in Appendix 1.

3.1 Leaf, stem and fruit symptoms

Phyllosticta ampelicida is most likely to be found on leaves, stems and fruit of grapevines and minor hosts. Sources could include imported fruit or cuttings.

In the vineyard, symptoms are visually most evident on leaves in the spring, leaves, stems and fruit in the summer, and on stems and fruit in autumn and winter.

Infected **leaves** can be easily detected by circular lesions which vary in size (Fig 2b) but are typically small under field conditions (Fig 2a and d), appear brown or tan with reddish margins, and characteristically contain small black pycnidia (Fig 2**Figure** c). When inspecting vines, leaf lesions may be observed associated with mummified berries from the previous or current seasons (Fig 2d).

On **shoots, petioles and tendrils,** the lesions are initially tan, brown turning purple to black, sunken, elliptical to elongated and contain pycnidia or pseudothecia observed as small black dots (Fig 3b and c). The bark may split, but stem infections remain localised, and do not usually extend more than several centimetres.




Figure 2. Symptoms of *Phyllosticta ampelicida* on leaves of *Vitis vinifera* (cv. Riesling); (a) on the vine, (b) on the laboratory bench, (c) under the dissecting microscope with pycnidia evident and (d) leaf lesions caused by conidia ejected from infected berries above. (Photos by M. Sosnowski, SARDI)




Figure 3. Symptoms of *Phyllosticta ampelicida* on stems and petioles of *Vitis vinifera* (cv. Riesling); (a) on young green shoots on the vine, (b) on the laboratory bench, (c) under the dissecting microscope, and (d) on older lignified shoots of a *Vitis* interspecific hybrid (cv. Vignoles) on the vine. (Photos a-c by M. Sosnowski, SARDI and photo d by W. Wilcox, Cornell University)

On the **developing fruit**, small, pale dots (approx. 1mm diam.) appear, rapidly becoming surrounded by a widening brown ring. If lesion expansion is halted due to an application of some fungicides or the development of age-related resistance in half grown berries, pale-brown coloured spots with a dark ring (bird's eye effect) and a sunken centre, about 6mm diam., may result (Fig 4a). More typically, however, a chocolate-brown lesion expands through the berry (Fig 4b and c) until it becomes completely rotted. Berries then shrivel and turn dark brown with numerous black pseudothecia or pycnidia developing over the surface (Fig 4d). Eventually, the fruit becomes dry and shrivelled, turning into hard, blue-black mummies that often remain firmly attached to the pedicel, although some may be shed from the vine (Fig 4e).




Figure 4. Symptoms of *Phyllosticta ampelicida* on grape berries of *V. labrusca* (cv. Concord); (a-c) progression from small lesions to mummified berries, (d) pseudothecia on surface of mummified berry under dissecting microscope and (e) mummified berries. (Photos a, b, d and e by M. Sosnowski, SARDI and c by W. Wilcox, Cornell University)

3.2 Confusion with other diseases

Symptoms can be confused with those of **black spot (anthracnose) disease** caused by the fungus *Elsinoë ampelina,* which is endemic in Australia. Leaf lesions are similar when viewed from a distance, but do not contain the small black pycnidia typically found in black rot lesions. Stems and petioles develop raised cankers with sunken centres, whereas black rot lesions are less "three-dimensional". Infected fruit develop similar lesions initially, but those of black spot do not usually expand to include the entire berry nor evolve into mummified fruit (Fig 5).




Figure 5. Symptoms of black spot (anthracnose) disease caused by the fungus *Elsinoë ampelina* on *Vitis vinifera* (table grape cv. Red globe); (a) leaf lesions, (b) stem cankers and (c) berry lesions. (Photos by M. Sosnowski, SARDI)

Symptoms of **Phomopsis cane and leaf spot**, which is caused by the fungus *Diaporthe ampelina*, may also be confused with those of black rot (Fig 6). Grapevines infected by *D. ampelina* develop leaf lesions that are smaller than those due to black rot; they are often surrounded with a translucent halo, and do not contain pycnidia. Stem lesions are sometimes smaller in size, but when larger, the two diseases are very difficult to distinguish on the basis of lesions on the green shoots. Lesions on lignified canes often have a bleached appearance and develop into larger basal cankers. Fruit symptoms can be similar to black rot, with shrivelled mummified berries; however, this phase of the disease does not usually occur with *D. ampelina* infections until shortly before harvest, whereas it typically occurs with veraison with black rot.

There is also a complex of pathogens, including *Alternaria* spp., *Aspergillus niger*, *Botrytis cinerea*, *Cladosporium* spp., *Penicillium* spp., and *Rhizopus arrhizus*, which may cause fruit rot of grapevines, and thus should also be considered when diagnosing black rot.



Figure 6. Symptoms of Phomopsis cane and leaf spot disease caused by the fungus *Diaporthe ampelina* on *Vitis vinifera* (cv. unknown); (a) small leaf lesions with halos, (b) green stem lesions, (c) bleached cane, (d) basal stem canker and (e) bunch symptoms. (Photos a-d by B. Rawnsley, SARDI and e by R. Emmett, DPI Victoria)

3.3 Sampling methods

As visual symptoms can be found on leaves, stems and fruit of grapevines, they should all be targeted during any inspections of vineyards, post-entry quarantine or examination of grapevine material at the border. In the vineyard, symptoms are likely to be most obvious on leaves first, but mummified fruit could also be present from infection in the previous season. Symptoms to look for are brown circular lesions (often containing small black pycnidia) with reddish margins on leaves; brown/purple elliptical or elongated lesions on stems; and brown circular lesions on fruit or shrivelled berries (details in section 3.1). In vineyards, symptoms can be present both during the growing season and dormancy, as infected, mummified fruit often remain hanging on the vine since they are not detached from the rachis as easily as healthy berries that merely shrivel on the vine through normal dehydration if not harvested at maturity.

Symptomatic leaf, stem and fruit samples should be sealed in a polyethylene bag and transported to a diagnostic facility, keeping the samples cool during the entire process. For national guidelines on response to an emergency plant pest refer to PLANTPLAN (PHA 2009).

4 IDENTIFICATION

Phllysticta ampelicida can only be confirmed by DNA sequence analysis, as morphology is unreliable for the identification of *Phyllosticta* species.

4.1 Isolation of pathogen from symptomatic plant material

Equipment

Dissecting microscope Scalpel Trays, polyethylene bags and paper towel Slides and coverslips Sterile distilled water (SDW)

Autoclave Sterile plastic petri dishes Laminar flow cabinet Incubator (25°C, white fluorescent lights) Scalpel and forceps Bunsen burner and alcohol Sodium hypochlorite solution (0.5%) Sterile distilled water (SDW) Sterile filter paper Parafilm

Media preparation:

Potato Dextrose Agar – half-strength (½ PDA) Hoffman et al. (2002)

Potato Dextrose Agar (PDA)	9.8 g
Bacto agar (granulated agarose)	3.5 g
Distilled water	500 ml

Mix ingredients and autoclave at 121°C for 15 minutes, mixing well before pouring into sterile plastic petri dishes.

Method:

Symptomatic leaf, stem and fruit material should be placed in a tray lined with wet paper towel and sealed in a plastic bag for incubation under high humidity. The samples can be left on the bench in normal laboratory conditions. Within 1-3 days, pycnidia on lesions may exude cirrhi or tendrils containing conidia which can be viewed under the dissecting microscope (Fig 7).




Figure 7. (a) Pycnidia of the sexual morph of *Phyllosticta ampelicida*, (b) exuding cirrhi on the leaf surface of grapevine cv. Riesling following incubation at high humidity overnight. (Photos by M. Sosnowski, SARDI)

In the laminar flow cabinet, remove sections of leaf, stem and fruit samples with lesions using a scalpel. Surface sterilise sections in 0.5% sodium hypochlorite for 60 s, rinse in sterile distilled water (SDW) and place on sterile filter paper to dry. Excise small pieces (approx. 3 mm x 3 mm) from lesions and place onto surface of agar media.

Alternatively, cirrhi exuding from pycnidia on lesions can be carefully removed with a flame-sterilised scalpel and placed directly onto agar surface.

Seal plates with parafilm strips and place into incubator at 25°C under continuous fluorescent light. Check plates for contamination at regular intervals and subculture putative colonies of *G. bidwellii* if necessary. Pure cultures should be incubated for 15-20 days in order to induce sporulation.

4.2 Morphological examination

Half-strength PDA is the optimal medium for identification and conidia production (Fig 8). Mycelium growth rate ranges from 1-2 mm/day and appears speckled and irregular. Pycnidia formation can occur within 7 days for some isolates and after 15 days all isolates produce fruiting bodies, from which conidia can be harvested.

Conidial ooze can be scraped from pycnidia on agar cultures or leaf lesions using a scalpel, and mounted on a glass slide in a drop of SDW with a cover slip placed on top.

Viewed under a compound light microscope, conidia appear as one celled, hyaline, broadly ovoid, ellipsoidal or almost globose, somewhat clavate when young and slightly indented; spores are $5-12 \times 4-7 \mu m$ in size (Fig 9).




Figure 8. *Phyllosticta ampelicida* cultures on ½ Potato Dextrose Agar; (a) 15 days and (b-f) 20 days of incubation at 25°C under continuous fluorescent light. (e & f) Under a dissecting microscope, conidia are observed oozing from pycnidia. (Photos by M. Sosnowski, SARDI)




Figure 9. (a) Pycnidia and (a & b) conidia of *Phyllosticta ampelicida* from cultures on ½ PDA incubated for 15 days at 25°C under continuous fluorescent light. (Photos by M. Sosnowski, SARDI)

4.3 Molecular method

Laboratory supplies

Protective gloves 2.0, 200, and 1000 μL sterile barrier pipette tips 2.0, 20, 200, and 1000 μL pipette tips Microcentrifuge Microcentrifuge tubes 1.5 mL 0.2 mL PCR tubes Thermocycler Gel tray with suitable comb/s, electrophoresis tank and powerpack UV transilluminator Camera/gel documentation system

4.3.1 DNA extraction

Genomic DNA (gDNA) may be extracted from pure fungal cultures and infected plant tissues using commercially available plant DNA extraction kits; e.g. QIAGEN Plant DNeasy kit, Bioline ISOLATE II Plant DNA extraction kit, etc. The preference is for gDNA to be eluted in nuclease-free H₂O. Store gDNA at -20°C.

4.3.2 Conventional PCR

Primers and expected product sizes

For amplification and sequencing of the ITS region:

V9G (forward)	5'-TTACGTCCCTGCCCTTTGTA-3' (de Hoog & Gerrits 1998)
ITS4 (reverse)	5'-TCCTCCGCTTATTGATATGC-3' (White et al. 1990)

Expected product size: approx. 750bp

For amplification and sequencing of partial region of the actin (ACT) gene (Carbone & Kohn 1999):

ACT512F (forward)	5'- ATGTGCAAGGCCGGTTTCGC -3'
ACT783R (reverse)	5'- TACGAGTGCCTTCTGGCCCAT -3'

Expected product size: approx. 280bp

For DNA sequencing analysis, it is preferred (but not essential) that a high-fidelity PCR master mix is used. Using the high-fidelity master mix will avoid the introduction of PCR errors into the product for DNA sequencing purposes, and to minimise pipetting errors in the preparation of PCR master mix. High-fidelity PCR master mixes can be purchased on a commercial basis, e.g. Phusion® High-Fidelity PCR Master Mix, Platinum® Taq DNA polymerase High Fidelity, etc.

DNA gel running buffer – TBE Buffer

Tris-Borate-EDTA (TBE) gel running buffer can be purchased commercially in concentrated liquid format. Follow the manufacturer's instructions to dilute it to a 1X concentration. Alternatively, it can be made up from the following components:

	Per 1L	Final conc.
Tris base	54.0 g	0.4 M
Boric acid	27.5 g	0.05 M
0.5M EDTA pH 8.0	20.0 mL	0.001 M

Dissolve components in 1L distilled water. Store at room temperature.

Note: Alternatively TAE (1×Tris-acetate-EDTA) buffer can be used.

Agarose gel

	Per 100 mL	Conc.
DNA grade agarose	1.0 g	1%
	1.5g	1.5%
TBE	100 mL	1x

Dissolve the molecular biology-grade agarose in TBE buffer in a heat-proof glass container (e.g. beaker or Schott bottle) by heating in a microwave. Once the bottle is slightly cool to the touch, pour into the gel tray with comb. It will take approximately 30 mins to set at room temperature (20-22°C).

WARNING: The container and the content are extremely hot. Handle the container with care using heatproof gloves.

DNA Loading dye and DNA stain

Certain dye powders are hazardous in concentrated forms. Therefore due to the workplace health and safety considerations of dye powders, DNA loading dye should be purchase from commercial companies.

Either Ethidium Bromide or GelRed® may be used to stain double-stranded DNA to visualise it on an agarose gel.

PCR

Prepare the following PCR mix as described below for the number of test samples, a blank, and one extra.

<u>Reagents</u>	<u>Volume (μL) per reaction</u>
Phusion ® Master Mix (2X)	12.5
Forward Primer (10µM)	0.5
Reverse Primer (10µM)	0.5
Nuclease-free H ₂ O	10.5
Total	24.0

Add 1 μ L of gDNA from pure fungal culture as template, place in a thermocycler using the following program:

	Initial denaturation	98°C	30 secs	1 cycle
--	----------------------	------	---------	---------

Denaturation	98°C	10 secs	30 cycles
Annealing	55 or 60°C*	30 secs	
Extension	72°C	30 secs	
Final Extension	72°C	5 mins	1 cycle
Cool	12°C	10 mins	1 cycle

* Anneal temperatures for the following PCR: ITS at 55°C and ACT at 60°C.

When the PCR is complete, mix 2μ l of each PCR sample with 1μ l of DNA loading dye.

Load the samples into a 1% TBE agarose gel.

Run the gel in TBE buffer at 100V for 40 minutes.

Visualise and photograph the gel under UV-light.

Sequence PCR product

Once it is confirmed that there is a single PCR product of expected size, prepare the PCR product for sequencing. Refer to the sequencing facility's guidelines for sample preparation and shipment.

DNA sequence analysis

Sequences from the ITS region and partial region of the ACT gene should only be compared to *P. ampelicida* ex-neotype strain ATCC 200578, GenBank accession KC193586 (ITS) and KC193581 (ACT). For a positive identification, the ITS sequence of the sample must be 100% match to the reference sequence, while ACT must be 99-100% match.

5 CONTACTS FOR FURTHER INFORMATION

Dr Mark Sosnowski

South Australian Research and Development Institute GPO Box 397 South Australia 5001 Ph: 08 8303 9489 Fax: 08 8303 9393 Email: mark.sosnowski@sa.gov.au

Prof Wayne Wilcox

Cornell University New York State Agricultural Experiment Station PO Box 462 Geneva, NY 14456 USA Ph: +1 315 787 2335 Email: wfw1@cornell.edu

6 ACKNOWLEDGEMENTS

This document is based on Information taken from PaDIL - Plant Biosecurity Toolbox which was originally authored by Gary Kong. Diagnostic methods for Black Rot of Grapes - *Guignardia bidwellii* [http://www.padil.gov.au/pbt] accessed 5 June 2009.

The information was reviewed and written into the V1 protocol by Dr Mark Sosnowski with funding from the Department of Agriculture and Water Resources (DAWR) through the Diagnostic Training Scholarship programme.

Wayne Wilcox and Judy Burr (Cornell University, NY, USA) assisted with training and providing details for the diagnostic protocol.

The molecular section of the protocol was developed by Yu Pei Tan, Qld DAF (2016).

The protocol was reviewed by Sri Kanthi de Alwis, Victorian DPI (2010) and Yu Pei Tan, Qld DAF (2016). The molecular section was reviewed by Jeyaseelan Baskarathevan, Plant Health and Environment Laboratory, MPI, New Zealand

7 REFERENCES

- Carbone, I., Kohn, L.M. (1999) A method for designing primer sets for speciation studies in filamentous ascomycetes. *Mycologia* 91(3): 553-556.
- Ferrin DM and Ramsdell DC (1977) Ascospore dispersal and infection of grapes by *Guignardia bidwellii*, causal agent of grape black rot disease. *Phytopathology*, 67: 1501-1505.
- Hoffman LE, Wilcox WF, Gadoury DA and Seem RC (2002) Influence of grape berry age on susceptibility to *Guignardia bidwellii* and its incubation period length. *Phytopathology*, 92: 1068-1076.
- de Hoog GS, Gerrits van den Ende AHG (1998) Molecular diagnostics of clinical strains of filamentous basidiomycetes. *Mycoses* 41:183–189
- PHA (2009) PLANTPLAN: Australian Emergency Plant Pest Response Plan. Version 1 May 2009, Plant Health Australia, Canberra ACT [http://www.planthealthaustralia.com.au/]
- Sivanesan A and Holliday P (1981) CMI Descriptions of Pathogenic Fungi and Bacteria, No. 710. Wallingford, UK: CAB International.
- Spotts RA 1977. Effect of leaf wetness duration and temperature on the infectivity of *Guignardia bidwellii* on grape leaves. *Phytopathology*, 67: 1378-1381.
- White TJ, Bruns T, Lee S, Taylor J (1990) Amplification and direct sequencing of fungal ribosomal RNA genes for phylogenetics. In: Innis MA, Gelfand DH, Sninsky JJ, White TJ (eds) PCR protocols: a guide to methods and applications. Academic Press, San Diego, pp 315–322
- Wilcox W (2003) Grapes: Black rot (*Guignardia bidwellii* (Ellis) Viala and Ravaz.) Cornell Cooperative Extension Disease Identification Sheet No. 102GFSG-D4, Cornell University.
- Zhang K, Zhang N, Cai L (2013) Typification and phylogenetic study of *Phyllosticta ampelicida* and *P. vaccinii*. Mycologia 105: 1030–1042.

7.1 Other useful references

- Alves A, Phillips AJL, Henriques I and Correia A (2007) Rapid differentiation of species of botryosphaeriaceae by PCR fingerprinting. *Research in Microbiology*, 158: 112-121.
- CAB International Distribution Maps of Plant Diseases, 1991, April (Edition 4), Map 81.
- Ferrin DM and Ramsdell DC (1978) Influence of conidia dispersal and environment on infection of grape by *Guignardia bidwellii*. *Phytopathology*, 68: 892-895.
- Hoffman LE and Wilcox WF (2002) Utilizing epidemiological investigations to optimize management of grape black rot. *Phytopathology*, 92: 676-680.

- Jailloux F (1992) *In vitro* production of the teleomorph of *Guignardia bidwellii*, causal agent of black rot of grapevine. *Canadian Journal of Botany*, 70: 254-257.
- Janex-Favre MC, Pargueyleduc A and Jailloux F (1993) The ontogeny of pycnidia of *Guignardia bidwellii* in culture. *Mycological Research*, 97: 1333-1339.
- Janex-Favre MC, Parguey-Leduc A and Jailloux F (1996) The ontogeny of perithecia in *Guignardia bidwellii*. *Mycological Research* 100: 875-880.
- Jermini M and Gessler C (1996) Epidemiology and control of grape black rot in southern Switzerland. *Plant Disease*, 80: 322-325.
- Kummuang N, Diehl SV, Smith BJ and Graves CH (1996) Muscadine grape berry rot diseases in Mississippi: Disease epidemiology and crop reduction. *Plant Disease*, 80: 244-247.
- Kummuang N, Smith BJ, Diehl SV and Graves CH (1996) Muscadine grape berry rot diseases in Mississippi: Disease identification and incidence. *Plant Disease*, 80: 238-243.
- Pezet R and Jermini M (1989) Black rot of grapevine: symptoms, epidemiology and control. *Revue Suisse de Viticulture, d'Arboriculture et d'Horticulture,* 21: 27-34.
- Ramsdell DC and Milholland RD, (1988) Black rot. In: APS, ed. Compendium of Grape Diseases. St Paul, USA: APS Press, 15-17.

Rao R and Kale SB (1965) A new species of *Guignardia* from India. *Mycopathologia*, 27: 13-14.

Spotts RA (1980) Infection of grape by *Guignardia bidwellii* - factors affecting lesion development, conidial dispersal, and conidial populations on leaves. *Phytopathology*, 70: 252-255.

8 APPENDICES

8.1 Disease cycle

The black rot disease cycle is illustrated in Fig 10. The black rot fungus overwinters primarily in mummified fruit within the vine and on the ground, although it also can overwinter for at least 2 years within lesions of infected shoots that are retained as canes or spurs. Spring rains trigger release of ascospores (airborne sexual spores) that form within mummies on the ground and in the trellis, and these can be blown for moderate distances by wind. Conidia (asexual spores) can also form, both within cane lesions or on mummies that have remained within the trellis, and these are dispersed short distances within the vine or to neighbouring vines by splashing rain drops. This process can start within 1 hour of the onset of rain (at least 0.3 mm) and may continue for up to 8 hours after rainfall ceases (Ferrin and Ramsdell 1977). Infection occurs when either spore type lands on susceptible green tissue such as leaves, blossoms and young fruit and it remains wet for 6-24 h depending on temperature (Spotts 1977).

Ascospores slowly germinate, often taking 36 to 48 hours, but eventually penetrate the young leaves and fruit stems (pedicels). The period of time required for symptoms to appear after the occurrence of infection depends on both the temperature and the age of the tissue at the time of infection, usually taking 8 to 25 days. The period during which these overwintering spores are available to cause infections depends on their source. From mummies on the ground, significant discharge of ascospores begins about 2 to 3 weeks after bud break and is virtually complete within 1 to 2 weeks after the start of flowering. In contrast, mummies within the trellis can continue to release both conidia and ascospores from the early pre-flowering period through veraison. From overwintering cane lesions, conidia can be dispersed from bud break through mid-summer.

Pycnidia, produced within the lesions, continue to release conidia during wet weather throughout the season. The conidia will germinate and infect leaves, blossoms and young fruit and can cause substantial spread of the disease under warm and rainy conditions, particularly if berries are still susceptible to infection after conidia develop. Most berries that become infected near the end of their period of susceptibility do not show symptoms until at least 3 weeks later, and the majority do not begin to rot until 4 to 5 weeks after the infection event. Fruit are most susceptible to infection by the fungus from mid-flowering to about 6 weeks after flowering, and become resistant to infection at maturity. Some older literature has reported that berries become resistant when they reach 5% to 8% sugar, while research in New York indicates that berries become resistant much earlier, 3 to 4 weeks after flowering.

In summary, spore production, dispersal, infection and continued disease are favoured by warm, humid conditions, summer rainfall and persistent dew. .




Figure 10. The disease cycle of black rot, caused by Phyllosticta ampelicida (Wilcox 2003)