

National Diagnostic Protocol

Uromyces viciae-fabae

The cause of lentil rust



NDP 31 V1

© Commonwealth of Australia

Ownership of intellectual property rights

Unless otherwise noted, copyright (and any other intellectual property rights, if any) in this publication is owned by the Commonwealth of Australia (referred to as the Commonwealth).

Creative Commons licence

All material in this publication is licensed under a Creative Commons Attribution 3.0 Australia Licence, save for content supplied by third parties, logos and the Commonwealth Coat of Arms.



Creative Commons Attribution 3.0 Australia Licence is a standard form licence agreement that allows you to copy, distribute, transmit and adapt this publication provided you attribute the work. A summary of the licence terms is available from <http://creativecommons.org/licenses/by/3.0/au/deed.en>. The full licence terms are available from <https://creativecommons.org/licenses/by/3.0/au/legalcode>.

This publication (and any material sourced from it) should be attributed as: Subcommittee on Plant Health Diagnostics (2015). National Diagnostic Protocol for *Uromyces viciae-fabae* – NDP31 V1 (Eds. Subcommittee on Plant Health Diagnostics) Authors Negussie, T, Pretorius, ZA, McTaggart, AR and Shivas, RG; Reviewer Toome, M. ISBN 978-0-9945112-1-8 CC BY 3.0.

Cataloguing data

Subcommittee on Plant Health Diagnostics (2015). National Diagnostic Protocol for *Uromyces viciae-fabae* – NDP32 V1. (Eds. Subcommittee on Plant Health Diagnostics) Authors Negussie, T, Pretorius, ZA, McTaggart, AR and Shivas, RG; Reviewer Toome, M. ISBN 978-0-9945112-1-8.

ISBN; 978-0-9945112-1-8

Internet

Report title is available at: <http://plantbiosecuritydiagnostics.net.au/resource-hub/priority-pest-diagnostic-resources/>

Department of Agriculture and Water Resources

Street Address 1: [18 Marcus Clarke Street](#), Canberra City ACT 2601

Street Address 2: [7 London Circuit](#), Canberra City ACT 2601

Postal Address: GPO Box 858, Canberra City ACT 2601

Switchboard Phone: 02 6272 3933

Web: <http://www.agriculture.gov.au>

Inquiries regarding the licence and any use of this document should be sent to: copyright@agriculture.gov.au.

The Australian Government acting through the Department of Agriculture and Water Resources has exercised due care and skill in the preparation and compilation of the information and data in this publication. Notwithstanding, the Department of Agriculture and Water Resources, its employees and advisers disclaim all liability, including liability for negligence, for any loss, damage, injury, expense or cost incurred by any person as a result of accessing, using or relying upon any of the information or data in this publication to the maximum extent permitted by law.

Purpose

National Diagnostic Protocols (NDPs) are diagnostic protocols for the unambiguous taxonomic identification of plant pests. NDPs:

- are a verified information resource for plant health diagnosticians
- are consistent with ISPM No. 27 – Diagnostic Protocols for Regulated Pests
- provide a nationally consistent approach to the identification of plant pests enabling transparency when comparing diagnostic results between laboratories; and,
- are endorsed by regulatory jurisdictions for use (either within their own facilities or when commissioning from others) in a pest incursion.

Where an International Plant Protection Convention (IPPC) diagnostic protocol exists it should be used in preference to NDPs although NDPs may contain additional information to aid diagnosis. IPPC protocols are available on the IPPC website:

<https://www.ippc.int/core-activities/standards-setting/ispms>

Process

NDPs are facilitated and endorsed by the Subcommittee on Plant Health Diagnostics (SPHD). SPHD reports to Plant Health Committee and is Australia's peak technical and policy forum for plant health diagnostics.

NDPs are developed and endorsed according to Reference Standards developed and maintained by SPHD. Current Reference Standards are available at

<http://plantbiosecuritydiagnostics.net.au/sphd/sphd-reference-standards/>

NDPs are living documents. They are updated every 5 years or before this time if required (i.e. when new techniques become available).

Document status

This version of the National Diagnostic Protocol (NDP) for *Uromyces viciae-fabae* is current as at the date contained in the version control box below.

PEST STATUS	Not present in Australia
PROTOCOL NUMBER	NDP 31
VERSION NUMBER	V1
PROTOCOL STATUS	Endorsed
ISSUE DATE	January 2016
REVIEW DATE	2021
ISSUED BY	SPHD

The most current version of this document is available from the SPHD website:

<http://plantbiosecuritydiagnostics.net.au/resource-hub/priority-pest-diagnostic-resources/>

Further information

Inquiries regarding technical matters relating to this project should be sent to:

sphds@agriculture.gov.au

Contents

1	INTRODUCTION.....	2
1.1	The Pathogen	2
1.2	Host range.....	2
1.3	Available protocols.....	2
2	TAXONOMIC INFORMATION.....	3
3	DETECTION	4
3.1	Plant parts affected	4
3.2	Symptom description	4
3.3	Diseases causing similar symptoms	6
4	IDENTIFICATION.....	7
4.1	Morphological methods.....	7
4.2	Molecular methods	9
5	CONTACTS FOR FURTHER INFORMATION.....	11
6	ACKNOWLEDGEMENTS.....	12
7	REFERENCES.....	13

1 INTRODUCTION

1.1 The Pathogen

Strains of the rust fungus *Uromyces viciae-fabae* show host specificity to legumes within the Fabaceae tribe of the Fabaceae (Barilli *et al.* 2011, Conner and Bernier 1982, Emeran *et al.* 2008, Hiratsuka 1933, Laundon and Waterston 1964). Systematic analysis indicated that strains of *U. viciae-fabae* either formed monophyletic groups (Emeran *et al.* 2008) or grouped according to their host species within one clade (Barilli *et al.* 2011). The rust causes partial defoliation and under heavy infection leads to premature plant death (Laundon and Waterston 1964, Negussie and Pretorius 2012). The rust fungus completes its full lifecycle of five spore stages on one host.

Uromyces viciae-fabae is present in Australia on faba beans, but has not been recorded on lentil.

1.2 Host range

The majority of investigations have reported that *U. viciae-fabae* has host specific strains on *Vicia* (faba bean), *Lathyrus* (sweet pea), *Lens* (lentil) and *Pisum* (pea). These included studies on host inoculations and morphology (Conner and Bernier 1982, Hiratsuka 1933) and molecular analyses (Barilli *et al.* 2011, Emeran *et al.* 2008). One study on the diversity of the rust in Japan determined that host specialization did not occur (Chung *et al.* 2004), although this work was based on a small number of isolates and has not been repeated.

1.3 Available protocols

Molecular diagnostic protocols have not been developed for *U. viciae-fabae* on lentil. The study by Chung *et al.* (2004) demonstrated that the Large Subunit (LSU) region of ribosomal DNA (rDNA) was not variable enough to distinguish between isolates from different hosts. Barilli *et al.* (2011) determined the ITS region differentiated lentil rust from other strains of *U. viciae-fabae*, however this sequence was not made available on GenBank.

2 TAXONOMIC INFORMATION

Kingdom:	Fungi
Phylum:	Basidiomycota
Class:	Pucciniomycetes
Order:	Pucciniales
Family:	Pucciniaceae
Genus:	<i>Uromyces</i>
Species:	<i>viciae-fabae</i> (Pers.) J. Schröt.
Synonyms:	<i>Aecidium leguminosarum</i> (Link) Rabenh. <i>Caeoma appendiculatum</i> Schltdl. <i>Caeoma leguminosarum</i> Link <i>Capitularia fabae</i> (Pers.) Syd. <i>Coemurus fabae</i> (Pers.) Kuntze <i>Nigredo fabae</i> (Pers.) Arthur <i>Puccinia fabae</i> Grev. <i>Puccinia fabae</i> (Alb. & Schwein.) Link <i>Puccinia fallens</i> Cooke <i>Puccinia globosa</i> Grev. <i>Puccinia polygona-avicularis</i> var. <i>fabae</i> Alb. & Schwein. <i>Trichobasis fabae</i> (Pers.) Lév. <i>Uredo fabae</i> Pers. <i>Uredo fabae</i> DC. <i>Uredo leguminosarum</i> Rabenh. <i>Uredo orobi</i> Schumach. <i>Uredo viciae-fabae</i> Pers. <i>Uromyces fabae</i> (Pers.) de Bary <i>Uromyces fabae</i> var. <i>orobi</i> (Schumach.) Jørst. <i>Uromyces orobi</i> (Schumach.) Lév. <i>Uromyces viciae</i> Fuckel

Common name: Lentil rust

3 DETECTION

3.1 Plant parts affected

All aerial plant parts are affected.

3.2 Symptom description

Negussie and Pretorius (2012) thoroughly described the symptoms from all spore stages. The aecial stage occurs on the abaxial surface of leaves and pods to form white aecial cups filled with orange-yellow spores. Uredinia then develop as dark brown pustules on both surfaces of the leaves, stems and pods (Figure 1-5). Finally telia are produced from the uredinia that are black in colour. A heavy infection will result in leaf drop and premature death.



Figure 1 & 2. Uredinia of *Uromyces viciae-fabae* on adaxial leaf surface and stems, Ethiopia. Images supplied by Tadesse Negussie and Zacharias Pretorius.



Figure 3. Uredinial infection by *Uromyces viciae-fabae* on adaxial leaf surface of partially resistant cultivar of lentil, Ethiopia. Image supplied by Tadesse Negussie and Zacharias Pretorius.



Figure 4. Uredinia of *Uromyces viciae-fabae* on adaxial leaf surface, Ethiopia. Image by Tadesse Negussie and Zacharias Pretorius.



Figure 5. Severe uredinal infection of stems and leaves by *Uromyces viciae-fabae*, Ethiopia. Image by Tadesse Negussie and Zacharias Pretorius.

3.3 Diseases causing similar symptoms

Uromyces viciae-fabae is the common cause of rust on lentil. Some species of rust on other host genera in the Fabaceae have indistinguishable symptoms from *U. viciae-fabae* on lentil, for example, strains of *U. viciae-fabae* on *Cicer arietinum* (chickpea), *Lathyrus* spp., *Lens* spp., *Pisum sativum* (pea), and *Vicia* spp. (Barilli *et al.* 2011). The leguminous genera *Lotus*, *Medicago*, *Pisum*, *Trifolium* and *Vigna*, are host to *Uromyces anthyllidis*, *U. striatus*, *U. pisi-sativi*, *U. trifolii-repentis* and *U. vignae*, respectively. These five species of *Uromyces* are similar in morphology to *U. viciae-fabae*, and all are present in Australia and many can be identified on the Rust Fungi of Australia Lucid Key (available: <http://collections.daff.qld.gov.au/web/key/rustfungi/Media/Html/browse.html>) (Shivas *et al.* 2014).

4 IDENTIFICATION

Identification of the lentil strain of *U. viciae-fabae* is based on host identity together with the morphology of teliospores and urediniospores. A specific molecular test has not been developed for the lentil strain of *U. viciae-fabae*, although sequences for many other strains of *U. viciae-fabae* are available for comparison on GenBank.

4.1 Morphological methods

The morphology of the lentil strain of *U. viciae-fabae* is considered identical to other strains of the pathogen. Identification of the host plant is important for the correct determination of the lentil strain. Common species of Fabaceae can be identified using this interactive key:

<https://gobotany.newenglandwild.org/dkey/fabaceae/>.

4.1.1 Microscopic identification

The urediniospores and teliospores can be removed from the leaf or stem surface with scalpel or forceps and mounted on a microscope slide in water or lactic acid. The slide should be heated and then examined with 100x oil immersion objective to visualise surface ornamentation and spore size.

Spermogonia mostly abaxial, amphigenous in small groups associated with aecia. **Aecia** mostly abaxial, predominantly along veins, surrounding spermogonia or scattered, peridium cupulate, white, 0.3–0.4 µm diam. **Aeciospores** 18–26 × 15–21 µm, broadly ellipsoid, hyaline, finely verruculose, wall 1–1.5 µm thick. **Uredinia** amphigenous, yellowish brown, 0.5 mm diam. **Urediniospores** are 22–32 × 17–25 µm, broadly ellipsoid, uniformly echinulate, with 3–5 germ pores equatorial or scattered, and a light golden brown wall 1–2.5 µm thick (Figure 6). **Telia** adaxial or amphigenous, exposed, blackish brown, compact, 1–2 mm diam. **Teliospores** are ellipsoidal, obovoidal or cylindrical, with a rounded or sub-acute apex; 24–40 × 17–26 µm; wall chestnut-brown, smooth, 1–3 µm thick at the sides and 5–12 µm at the apex; pedicels brownish, up to 100 µm long (Figure 7) (from Negussie and Pretorius 2012).




Figure 6. Urediniospores and teliospores of *Uromyces viciae-fabae* on *Vicia sativa*, Tasmania (BRIP 60149). Image by Alistair McTaggart.



Figure 7. Teliospores of *U. viciae-fabae* on *Vicia sativa*, Tasmania (BRIP 60149). Image by Alistair McTaggart.

4.2 Molecular methods

4.2.1 Molecular barcoding of *U. viciae-fabae*

Amplified copies of the Large Subunit (LSU) region of rDNA can be sequenced and compared to known sequences on GenBank for identification of rust fungi. The LSU region is more easily sequenced than the ITS region for rust fungi, as the ITS may contain indels that inhibit direct sequencing. The ITS2-LSU region can be amplified with primers *Rust 2INV* and *LR7*. In cases when a product is not amplified or is of low concentration, a nested reaction can be performed using the primers *LROR* and *LR6*. In the case of *U. viciae-fabae*, sequences of the ITS region are unavailable for comparison with the lentil strain.

Equipment and reagents

- Thermocycler
- Taq polymerase and PCR components
- Micropipettes and aerosol resistant tips
- Disposable gloves (powder free)
- Gel electrophoresis apparatus

Primers (LSU primers):

Rust 2INV: 5'- GATGAAGAACACAGTGAAA -3' (Aime, 2006)

LR7: 5'- TACTACCACCAAGATCT -3' (Vilgalys and Hester, 1990)

LROR: 5'- ACCCGCTGAACTTAAGC -3' (Vilgalys and Hester, 1990)

LR6: 5'- CGCCAGTTCTGCTTACC -3' (Vilgalys and Hester, 1990)

DNA Extraction

Any standard fungal DNA extraction protocol can be used for rust fungi. A recommended protocol for DNA extraction from rust fungi is the UltraClean Microbial DNA Isolation Kit (MoBio Laboratories, Solana Beach, CA, USA)..

Between 5 and 20 rust sori are excised from a leaf using fine forceps or a scalpel and placed into extraction buffer. The kit protocol is then followed to completion and DNA is stored at -20°C.

DNA amplification protocol

1. Prepare PCR cocktail on ice in a sterile microcentrifuge tube. This reaction can be performed with any polymerase enzyme according to the manufacturer conditions.

For each 25 µl sample the cocktail will contain:

PCR buffer (10x)	2.500 µL (final concentration 1x)
MgCl ₂ (50 Mm)	0.750 µL (final concentration 1.5mM)
dNTPs (10 mM)	0.500 µL (final concentration 200µM)
Forward primer (10 µM)	0.250 µL (final concentration 0.2µM)
Reverse primer (10 µM)	0.250 µL (final concentration 0.2µM)
Taq	0.200 µL (final concentration 1%)
H ₂ O	20.550 µL

To prepare the cocktail, multiply the above volumes by the number of samples and add to a single tube, 24 μ L aliquots are then made into 0.2 mL tubes.

2. Add 1 μ L of DNA template to 24 μ L of PCR cocktail
3. Run PCR

Cycle conditions:

First reaction with primers *Rust 2INV* and *LR7*

Denaturation	94°C for 4 minutes	x 1 cycle
Denaturation	94°C for 30 sec	} x 45 cycles
Annealing	57°C for 45 sec	
Extension	72°C for 1.5 min	
Final extension	72°C for 7 min	x 1 cycle

4. Run a 5 μ L aliquot of this reaction on a 1–1.5% agarose gel to confirm successful amplification. A ~1200 base pair product should be expected for the reaction with *Rust 2INV* and *LR7*. If a product is observed proceed to step 6. If weak or no product is observed continue with the following reaction:

Dilute 1 μ L of PCR product in 99 μ L of sterile H₂O. 1 μ L of this dilution is then used as template for the next reaction with primers *LROR* and *LR6* with the protocol from step 1.

Nested reaction with *LROR* and *LR6*

Denaturation	94°C for 2 min	x 1 cycle
Denaturation	94°C for 30 sec	} x 45 cycles
Annealing	59°C for 30 sec	
Extension	72°C for 1.5 min	
Final extension	72°C for 7 min	x 1 cycle

5. Run a 5 μ L aliquot of the nested reaction on a 1–1.5% agarose gel to confirm successful amplification. A ~1000 base pair product should be expected for the nested reaction with *LROR* and *LR6*.
6. Successful PCR product should be sequenced. An example of a third party sequencing company is Macrogen, Korea. The directions for sample submission of the third party should be followed.
7. Sequences should be determined using chromatograms from both primers. A comparison of the sequence should be made with sequences of *U. viciae-fabae* on GenBank using a nucleotide BLAST search (<http://blast.ncbi.nlm.nih.gov/Blast.cgi>). A 99-100% sequence identity to any of the following *U. viciae-fabae* LSU sequences AB115592-AB115611, AY745695, KJ716343 indicates the specimen is *U. viciae-fabae* and the host must be identified to confirm it is the lentil strain. High sequence identity to HQ317516 *U. phaesoli*, which likely is a misapplication of this name on *Pisum*, should also be considered a positive identification of *U. viciae-fabae*.

5 CONTACTS FOR FURTHER INFORMATION

Zacharias A. Pretorius
Professor (Plant Pathology): Department of Plant Sciences
PO Box 339, Bloemfontein 9300, Republic of South Africa
T: +27 (0)51 4012466
E: PretorZA@ufs.ac.za

Dr Alistair McTaggart
E: alistair.mctaggart@gmail.com

Dr Roger Shivas
Biosecurity Queensland
Department of Agriculture and Forestry
DAF Level 2C East, Ecosciences Precinct, Basement 3 Loading Dock off Joe Baker Street,
Dutton Park, Qld 4102
T: 61 7 3255 4378
E: roger.shivas@daf.qld.gov.au

6 ACKNOWLEDGEMENTS

This protocol was developed by Drs Tadesse Negussie, Zacharias Pretorius, Alistair McTaggart and Roger Shivas.

The protocol was reviewed and verified by Dr Merje Toome, MPI, New Zealand.

7 REFERENCES

- Aime, M. C. 2006. Toward resolving family-level relationships in rust fungi (Uredinales). *Mycoscience*, 47, 112-122.
- Barilli, E., Satovic, Z., Sillero, J. C., Rubiales, D. & Torres, A. M. 2011. Phylogenetic Analysis of *Uromyces* Species Infecting Grain and Forage Legumes by Sequence analysis of Nuclear Ribosomal Internal Transcribed Spacer Region. *Journal of Phytopathology*, 159, 137-145.
- Chung, W. H., Tsukiboshi, T., Ono, Y. & Kakishima, M. 2004. Phylogenetic analyses of *Uromyces viciae-fabae* and its varieties on *Vicia*, *Lathyrus*, and *Pisum* in Japan. *Mycoscience*, 45, 1-8.
- Conner, R. L. & Bernier, C. C. 1982. Host range of *Uromyces viciae-fabae*. *Phytopathology*, 72, 687-689.
- Emeran, A. A., Roman, B., Sillero, J. C., Satovic, Z. & Rubiales, D. 2008. Genetic variation among and within *Uromyces* species infecting legumes. *Journal of Phytopathology*, 156, 419-424.
- Hiratsuka, Y. 1933. Studies on *Uromyces fabae* and its related species. *Japanese Journal of Botany*, 6, 329-379.
- Laundon, G. F. & Waterston, J. M. 1964. *Uromyces viciae-fabae*. *CMI Descriptions of Pathogenic Fungi and Bacteria*. CAB International.
- Negussie, T. & Pretorius, Z. A. 2012. Lentil rust: Present status and future prospects. *Crop Protection*, 32, 119-128.
- Shivas, R. G., Beasley, D. R. & McTaggart, A. R. 2014. Online identification guides for Australian smut fungi (Ustilaginomycotina) and rust fungi (Pucciniales). *IMA Fungus*, 5, 195-202.
- Vilgalys, R. & Hester, M. 1990. Rapid genetic identification and mapping of enzymatically amplified ribosomal DNA from several *Cryptococcus* species. *Journal of Bacteriology*, 172, 4238-4246.